

APPLICATION NOTE

Pharmaceuticals – TG 209 F1 Libra® -GC-MS

TG-GC-MS Measurements on Acetylsalicylic Acid (ASA)

Sonja Eichholz

Measurement Parameters


TG 209 F1 Libra® - GC-MS

TGA Parameters:

RT to 500°C
10 K/min
50 ml/min N₂
Sample mass: 1.71 mg
Al₂O₃ crucibles

GC Parameters:

Quasi-continuous
GC furnace: 250°C
Split: 10:1
Column: HP 5 ms, 30 m
Valve switching: 30 sec



1 TG 209 F1 Libra® with GC-MS coupling


Introduction

ASA is one of the most widely used non-opioid analgesics. By esterification of the phenolic hydroxyl group of salicylic acid with acetic acid, not only is a better local tolerance achieved, but also a stronger antipyretic, antiphlogistic and especially platelet aggregation inhibiting effect [1].

The common name is 2-acetoxybenzoic acid. It forms white, needle-shaped crystals which smell faintly of acetic acid. ASA is produced by means of the Kolbe-Schmitt synthesis, an electrophilic aromatic substitution [2].


[1] Mutschler, Arzneimittelwirkungen

[2] Laue, Reaktionsmechanismen



2 Chemical structure of acetylsalicylic acid

APPLICATIONNOTE TG-GC-MS Measurements on Acetylsalicylic Acid (ASA)



3 TGA results (black), DTA results (green), c-DTA[®] signal (red) and total ion chromatogram of the ASA measurements (purple)


Measurement Results

Simultaneous coupling of thermogravimetric analysis (TGA) and GC-MS (gas chromatography and mass spectrometry) enables easy analysis of the thermal behavior and decomposition products.

The melting of ASA can be detected by means of c-DTA[®] at 142°C. ASA decomposes in two steps (DTG minima at 175°C and 350°C). The total ion chromatogram reflects the

completely measured spectra as a function of temperature. Corresponding to DTG, two maxima can clearly be seen.

The 1st TG step can be attributed to the formation of acetic acid (mass numbers 43 and 60 figure 4); the 2nd step can be attributed to phenol (mass number 94) as a decomposition product.



4 TGA results (black) and total ion chromatogram of mass numbers 43, 60 and 94