

Industrial use of glycerol and the need for tribological model system testing



Industrial applications of Glycerol

- Glycerol is used in pharmaceutical, food and personal care industry^{1,2}.
- Typical applications are:
 - Used as a sweetener¹
 - Used to modify food viscosity¹
 - Used as a lubricant in oral care and personal care products²
 - Contact lens solutions for hydration
 - Biocompatible gels for transdermal patches

Tribological testing

- In this study, a tribological model system testing approach is introduced.
- It is used to compare tribological system behaviour with the rheological properties of different glycerol-water solutions..

Fig. 1 Tribological model system testing is carried out to study the lubricating properties of glycerol in personal care products.

Glycerol-water solutions

- Water-glycerol solutions with different glycerol concentrations were tested.
- The glycerol-water mixing ratios were:

H ₂ O : Glycerol mixing ratio
1 : 0
0.75 : 0.25
0.5 : 0.5
0.25 : 0.75
0 : 1

Tab. 1 Water-glycerol solutions mixing ratios.

Tribological model system testing

- A Kinexus pro+ rotational rheometer equipped with a Tribology cell was used in combination with a round-on-flat geometry.
- The balls were made from stainless steel, measuring 12.7 mm in diameter.
- The discs were made from a silicone elastomer and glued to the sample holder using an adhesive.
- Measurements were carried out in the form of Stribeck curves.
- New specimen were used for each test.
- Measurements were carried out at 20 °C.




Fig.2 Schematic of the concentric cylinder geometry.

Complementary shear viscosity measurements

- A Kinexus pro+ rotation rheometer coupled with cone-plate geometries with 60 mm and 65 mm respectively in diameter and 1° angle.
- Measurements were carried out at 20 °C.
- Shear rates in the range from 10 to 10³ s⁻¹ were covered.




Fig. 3 Schematic of the cone-plate geometry.

Differences in the tribological behavior due to shear viscosity of the glycerol-water solutions




Fig. 4 Stribeck curves.




Fig. 5 Viscosity-corrected Stribeck curves.

Stribeck curve measurements

- Coefficient of friction as a function of sliding speed and viscosity-corrected sliding speed (i.e. Stribeck curves) are shown for multiple glycerol-water ratios (see Fig. 4).
- The curves approximately collapse onto a master Stribeck curve (see Fig. 5).
- The results from this study suggest that the differences in the tribological behaviour found between the water-glycerol solutions can mainly be explained by differences in shear viscosity.

Shear viscosity measurements

- The glycerol solutions showed Newtonian flow behaviour.
- Steady-state shear viscosity data from the entire shear rate range were used to calculate the average data below (see Tab. 2).

H ₂ O : Glycerol mixing ratio	Average shear viscosity [Pas]	Average shear viscosity [Pas]
1 : 0	0.0013	0.0004
0.75 : 0.25	0.0021	0.0009
0.5 : 0.5	0.0064	0.0012
0.25 : 0.75	0.0230	0.0028
0 : 1	0.8259	0.0392

Tab. 2 Average steady-state shear viscosity of different water-glycerol solutions

Conclusion and outlook

- It was possible to distinguish between different Newtonian solutions with a reasonable degree of accuracy.
- Higher glycerol content provided lower coefficient of friction values at lower sliding speeds.
- These results show the importance of formulation in the food and personal care industry, where tribological characterization can help to study product perception.

References

- Chilakamarry, C.R. et al. (2021). Glycerol waste to value added products and its potential applications. Systems Microbiology and Biomanufacturing, 1, 378-396.
- Goyal, S., Hernández, N.B. and Cochran, E.W. (2021). An update on the future prospects of glycerol polymers. Polymer International, 70, 911-917.